PROTEÍNAS TOTALES

MÉTODO BIURET

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Gornall, A.G., Bardawill, C.J., David, M.M. (1948). Biol. Chem., 177, 751-766.

Peters, T. (1968). Clin. Chem., 14, 1147-1159.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

TOTAL PROTEIN

BIURET METHOD

For "in vitro" determination of total protein in serum or plasma

PRINCIPLE

In alkaline solution, the proteins form with copper(II) ions a colored complex, highly stable, which is spectrophotometrically measurable and proportional to the concentration of protein in the sample.

DIAGNOSTIC USE

The determination of total protein in serum is used in the study of the nutritional status or edematous processes.

Total protein levels are high (> 9.0g/dL) in: Hyperimmunoglobulinaemia, mono or polyclonal gammopathy, dehydration, chronic liver disease and and neoplasms, especially myeloma

Values lower than 6.0 g/dL are associated with protein loss in cases of gastroenteropathies, burns, nephrotic syndrome or due to decreased protein synthesis in chronic liver disease, malabsorption syndrome, agammaglobulinemia or malnutrition.

In cases of pregnancy, administration of intravenous fluids, chronic alcoholism, heart failure or hyperthyroidism, the blood protein values are also lower than normal.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 3 x 100 mL. (Ref. 99 71 80). Contents:

- A. 3 x 100 mL Biuret reagent
- B. 1 x 5 mL Standard

Ref. 99 96 02
Ref. 99 06 85

WORKING REAGENT PREPARATION

Reagent and standard are ready to use

REAGENT COMPOSITION

The reagent composition is as follows:

| | |
|-------------------------------|---------|
| NaOH | 0.47 M |
| Potassium iodide | 23.3 mM |
| Copper (II) sulphate | 6.5 mM |
| Sodium-potassium tartrate | 22.1 mM |
| Preservatives and stabilizers | |

Standard: Aqueous solution of protein equivalent to 5 g/dL (50 g/L).

STORAGE AND STABILITY

When kept at room temperature ($\leq 25^{\circ}\text{C}$), the reagent will remain stable until the expiration date stated on the label.

The Standard, however, should be stored at $2^{\circ}\text{-}8^{\circ}\text{C}$.

It is therefore recommended, to take the standard out of the kit box and store it accordingly.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent.
Working reagent blank >0.300

ADDITIONAL EQUIPMENT

General laboratory equipment.
Spectrophotometer or photometer with a 540 nm filter. Cuvette: 1 cm light path

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.
The safety statements are on the label. It is advisable to look at the MSDS before using the reagent.
The disposal of the residues has to be made according to legal local regulations.

SAMPLE

Serum or plasma, free from haemolysis. The sample will remain stable for up to 5 days, when kept at $2\text{-}8^{\circ}\text{C}$.
When using other corporal fluids as sample (amniotic liquid, urine, exudates, etc) it is advisable to take into account the margin values to test and the sensitivity of the reagent. See Performance characteristics.

PROCEDURE

| Technique | BL | ST | SA |
|-----------|------|------|------|
| | mL | mL | mL |
| Sample | --- | --- | 0.02 |
| Standard | --- | 0.02 | --- |
| Reagent | 1.00 | 1.00 | 1.00 |

Mix well and let stand for 10 min at room temperature ($20\text{-}25^{\circ}\text{C}$).

Reading

Wavelength: 540 nm
Blank: BL contents
Colour stability: a minimum of 3 hours

CALCULATIONS

$$\frac{\text{SA Abs}}{\text{ST Abs}} \times 5 = \text{g of protein / dL}$$

Where:

SA Abs: Sample absorbance
ST Abs: Standard absorbance

SI Units

$$(\text{g/dL}) \times 10 = \text{g/L}$$

REFERENCES VALUES

Adults: 6.4-8.3 g/dL
Children
Newborns: 4.6-7.0 g/dL
< 1 year: 5.1 - 7.3 g/dL
1 - 2 years: 5.6 - 7.5 g/dL
> 3 years: 6.0 - 8.0 g/dL

The stated values are for guidance. Each particular laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Performance of the reagent depends on the reagent itself and also depends on method and analyzer used.
The results indicated are obtained using a manual method.

Sensitivity, as detection limit: 0.10 g/dL
Linearity: Up to 12g/dL. For higher values, it is recommended to dilute the sample with saline (NaCl 0.9%) and assay once again. Multiply the final result by dilution factor.
Accuracy: 98.7 %.
Repetitivity as Variation Coefficient: 0.85%
Reproducibility, as Variation Coefficient: 1.13%
Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.
Details of the performance studies are available on request

INTERFERENCES

No interferences are known
The use of disposable equipment is recommended to prevent unwanted contamination.

QUALITY CONTROL

It is recommended to include control sera, Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85), in each test series for assuring the obtained values .

Each particular laboratory should establish its own quality control program and measures for avoiding results deviations

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

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Gornall, A.G., Bardawill, C.J., David, M.M. (1948). Biol. Chem., 177, 751-766.

Peters, T. (1968). Clin. Chem., 14, 1147-1159.

Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

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GB - TOTAL PROTEINS (BIURET), REAGENT

Warning

Hazard: Causes skin irritation. Causes serious eye irritation. Harmful to aquatic life with long lasting effects.

Precautionary: Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose the contents/container according to point 13 of the Safety Data Sheet.

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