
### PRINCIPIO

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# UREA U.V.


## UREASE – GLDH METHOD

For "in vitro" determination of urea in serum, plasma or urine

### PRINCIPLE

The hydrolysis of the urea present in the sample is catalyzed by urease producing ammonium and carbonate ions. Ammonium ions react with  $\alpha$ -ketoglutarate by the action of glutamate dehydrogenase, oxidizing the NADH to NAD<sup>+</sup>

The urea concentration present in the sample is proportional to the decrease in the concentration of NADH in the reaction.

### DIAGNOSTIC USE

Urea is the end product of protein metabolism; it is produced in the liver and eliminated from the blood in the kidneys. The serum level of urea is used in conjunction with creatinine, for the evaluation of renal function.

Usually the elevated levels of blood urea reflect some alteration of the excretory function of the kidney. This increase may also be due to poor liver function, cardiac failure or diet itself. Dehydration and bleeding can raise levels of blood urea.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

Kit 12 x 16 mL. (Ref. 99 13 05). Contents:

A. 12 vials of Freeze-dried enzymes	Ref. 99 16 00
B. 2 x 100 mL Buffer solution	Ref. 99 27 11
C. 1 x 5 mL Standard	Ref. 99 02 41

### WORKING REAGENT PREPARATION

Dissolve the contents of the enzymes vial A with the volume of buffer B stated on the label. Mix gently.

Standard is ready-to-use.

### WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Tris-HCl buffer pH 7.8	100 mM
Sodium $\alpha$ -ketoglutarate	6 mM
ADP	2 mM
NADH	0.20 mM
EDTA	4 mM
Urease	$\geq 4,000$ U/L
GLDH	$\geq 9,000$ U/L

Stabilizers and preservatives

Standard: Aqueous solution of Urea equivalent to 40 mg/dL. (6.6 mmol/L).

### STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The rehydrated reagent is stable up to 6 weeks at 2-8°C and 10 days at room temperature ( $\leq 25^\circ\text{C}$ ).

### Signs of reagent deterioration:

Presence of particles or turbidity. Working reagent blank < 1.00

### ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C. Cuvette: 1cm light-path

### SAMPLE

Serum, plasma or urine. Urea will remain stable in serum for at least 1 day at room temperature ( $\leq 25^\circ\text{C}$ ), 4-5 days at 2-8°C and 6 months when frozen ( $-20^\circ\text{C}$ ). In urine, urea will remain stable, 4-5 days when kept at 2-8°C, provided that the pH value be lower than 4.

If a urine sample is to be assayed, it should be previously diluted 1/100 with deionized water. Multiply the final result by 100.

### CAUTION

The reagents contain Sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the MSDS before using the reagent.

The disposal of the residues has to be made according to legal local regulations.

### QUALITY CONTROL

It is recommended to include control sera, Seriscann Normal (Normal Control Serum, Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum, Ref. 99 46 85), in each test series for assuring the obtained values.

### AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

### PROCEDURE

Let stand the reagents for 2-3 min. at 37°C.

Technique	ST mL	SA mL
Sample	--	0.01
Standard	0.01	--
Working reagent	1.00	1.00

Mix well and transfer to the measuring cuvette. Read the absorbance at 30s (Abs1) and after 1 min, read again (Abs2).

### Reading

Wavelength: 340 nm

Blank: Deionized water

Cuvette: 1 cm light-path

### CALCULATIONS

Determine  $\Delta\text{Abs}$  for each sample and the standard:

$$\Delta\text{Abs} = \text{Abs1} - \text{Abs2.}$$

$$\frac{\Delta\text{Abs SA}}{\Delta\text{Abs ST}} \times 40 = \text{mg urea / dL}$$

Where:

$\Delta\text{Abs SA}$ : Absorbance increase of sample

$\Delta\text{Abs ST}$ : Absorbance increase of standard

### S.I. Units

(mg/dL Urea) x 0.1665 = mmol/L Urea

### REFERENCES VALUES

Serum: 10 - 50 mg/dl Urea; 1.7 - 8.3 mmol/L Urea

Urine: 20 - 35 g Urea/24 h.; 333 - 583 mmol Urea /24 h

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

### Results as BUN (Blood Urea Nitrogen)

Mg/dL Urea x 0.467 = mg/dL BUN

### PERFORMANCE CHARACTERISTICS

Performance of the reagent depends on the reagent itself and also depends on method and analyzer used. The results indicated are obtained using a manual method.

Sensitivity, as detection limit: 2.0 mg/dL

Linearity: Up to 300 mg/dL of Urea.

Accuracy: 98.6 %.

Repetitivity as Variation Coefficient: 1.46%

Reproducibility, as Variation Coefficient: 1.65%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request

### INTERFERENCES

Fluoride as well as ammonium heparinate inhibit the reaction. Serum samples should be free from hemolysis and turbidity.

The use of disposable labware is recommended to prevent unwanted contamination.

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