LIPASA

MÉTODO COLORIMÉTRICO

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Lott, J. A. et al., (1986). Clin. Chem., 32, 1290 - 1302.

Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

LIPASE

COLORIMETRIC METHOD

For the "in vitro" determination of lipase in serum or plasma

PRINCIPLE

At alkaline pH the lipase substrate 1,2-O-dilauryl-rac-glycero-3- glutaric acid-(6-methylresorufin)-ester is cleaved by the catalytic action of pancreatic lipase to form 1,2-O-dilauryl-rac-glycerol and unstable compound glutaric acid-(6-methylresorufin)-ester. In that alkaline solution, this compound degrades to glutaric acid and methyl-resorufin. The color intensity of the red dye formed is directly proportional to the lipase activity. This activity can be measured at 578nm.

DIAGNOSTIC USE

The lipase is increased in pancreatitis from any source, and these values remain high longer than amylase values.

In chronic disorders of the pancreas the Lipase value is high and the same applies in cases of primary biliary cirrhosis, chronic renal failure and hemodialysis patients.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 80 mL. (Ref. 99 11 15). Contains:

A. 1 x 50 mL Buffer solution

Ref. 99 11 20

B. 1 x 30 mL Substrate

Ref. 99 12 15

C. 1 x 1mL Freeze-dried Standard

Ref. 99 13 07

The concentration is stated on the label of the vial.

WORKING REAGENT PREPARATION

Reagents are ready to use.

Rehydrate the standard with 1 mL of deionised water.

REAGENT COMPOSITION

The concentrations in the reagents solutions are:

A. Buffer solution

Bicin buffer pH 8.0

50 mM

Colipase

≥ 1 mg/L

sodium deoxycholate

1.8 mM

calcium chloride

12 mM

B. Substrate

Tartrate buffer pH 4.0

12 mM

1,2-O-dilauryl-rac-glycero-3-glutaric acid-

-(6-methylresorufin) - ester

0.27 mM

Taurodeoxy-cholate

9.0 mM

Surfactants and preservatives

STORAGE AND STABILITY

When stored at 2°-8°C, the components of the kit will remain stable until the expiration date stated on the label.

The standard, once rehydrated, is stable for 15 days at 2-8°C.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≥0.500.

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer, automated analyzer or photometer with thermostated cuvette.

SAMPLE

Serum or heparinized plasma. Use samples free from hemolysis. Do not use plasma with EDTA.

The lipase activity in serum is stable for 4 days at 2-8°C and for 24h at room temperature (≤ 25°C).

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the MSDS before using the reagent.

The disposal of the residues has to be made according to local regulations.

PROCEDURE

Bring reagents and the analyzer to 37°C.

Technique	BL mL	SA mL	ST mL
Water	0.01	---	---
Sample	---	0.01	---
Standard	---	---	0.01
Reagent A	1.00	1.00	1.00
Mix and incubate 5min, at 37°C. Then add Reagent B			
Reagent B	0.60	0.60	0.60

Mix well. After 60s read (Abs1). After 90s more, read again (Abs2)

Reading

Wavelength: 578nm

Blank: deionised water

Cuvette thermostated: 1cm light-path

CALCULATION

Determine ΔAbs for each sample and the standard:

ΔAbs = Abs2 - Abs1

$$\frac{\Delta\text{Abs SA} - \Delta\text{Abs BL}}{\Delta\text{Abs ST} - \Delta\text{Abs BL}} \times \text{Conc. ST (U/L)} = \text{U/L}$$

Results as μkat/L

U/L x 0.0167 = μkat/L

REFERENCE VALUES

Serum, plasma: < 38 U/L.

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method.

Sensitivity, as detection limit: 6.9 U/L

Linearity: Up to 300 U/L. For higher values, it is recommended to dilute the sample 1/5 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 5.

Accuracy: 104%

Within-run Precision as Coefficient of Variation: 3.01%

Run-to-run Precision as Coefficient of Variation: 3.65%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request

INTERFERENCES

Hemolysis interferes with the assay.

Concentrations of hemoglobin higher than 125 mg/dL, of bilirubin higher than 20 mg/dL and triglycerides higher than 625 mg/dL can interfere with the assay.

To perform this test, the use of disposable labware is highly recommended.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

AUTOANALYZERS

QCA has adaptation to other autoanalyzers on request.

REFERENCES

Rick, W., (1969), Zeitschrift Clin. Chem. Clin. Biochem., 7, 530 - 536.

Ziegenhorn, J. et al., (1979), Clin. Chem., 25, 1067.

Neumann, U., Kaspar, P., Ziegenhorn, J., (1984) Meth. Enz. Anal (3rd edition), 26 - 34.

Lott, J. A. et al., (1986). Clin. Chem., 32, 1290 - 1302.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

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