# COLESTEROL LÍQUIDO

## MÉTODO CHOD – POD

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Zoppi, F., Fellini, D. (1976), Clin. Chem., 22, 690-691.

Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

Zapico E. y Ordóñez J. (2002), Clin Invest Ar ter ioscl; 14(5):272-6.


# CHOLESTEROL LIQUID

CHOD - POD METHOD

For "in vitro" determination of cholesterol in serum or plasma.

## PRINCIPLE

Cholesterol present in serum or plasma, originates a coloured complex, according to the following reactions, which can be quantified by spectrophotometry.

## DIAGNOSTIC USE

The cholesterol in serum or plasma can be of exogenous origin, ingested with the diet, or endogenous synthesized primarily in the liver. It is transported by lipoproteins and it is excreted in the bile.

The study of the serum cholesterol level allows the detection and classification of the hyperlipidemia and risk assessment for heart disease.

Cholesterol levels are low in hypolipoproteinemia, hyperthyroidism and in some anemia.

Single test result cannot be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

## REAGENTS

**Kit 1 x 100 mL. (Ref. 99 52 82).** Contents:

A. 1 x 100 mL Reagent Ref. 99 52 20  
B. 1 x 5 mL Standard Ref. 99 02 48

**Kit 3 x 100 mL. (Ref. 99 52 80).** Contents:

A. 3 x 100 mL Reagent Ref. 99 52 20  
B. 1 x 5 mL Standard Ref. 99 02 48

**Kit 2 x 250 mL. (Ref. 99 50 12).** Contents:

A. 2 x 250 mL Reagent Ref. 99 01 59  
B. 1 x 5 mL Standard Ref. 99 02 48

## WORKING REAGENT PREPARATION

Reagent and standard are ready to use.

## WORKING REAGENT COMPOSITION

Concentration in the reagent solution is:

Mes buffer pH 6.5	75 mM
Phenol	6 mM
2,4-dichlorophenol	0.2 mM
4-aminoantipyrine	0.5 mM
Cholesterol esterase	≥ 500 kU/L
Cholesterol oxidase	≥ 300 kU/L
Peroxidase	≥ 1200 kU/L
Non reactive Stabilizers	

**Standard:** Solution of cholesterol in isopropanol/water equivalent to 200 mg/dL (5.18 mmol/L).

## STORAGE AND STABILITY

The components of the kit, when stored at 2-8°C, will remain stable until the expiration date stated on the label.

## Signs of reagent deterioration:

Presence of particles or turbidity in the reagent.  
Working reagent blank >0.500

## ADDITIONAL EQUIPMENT

General laboratory equipment.  
Spectrophotometer or photometer thermostable at 37°C.

## SAMPLE

Serum, or plasma. Samples are stable 1 week at 2-8°C, and up to three months at -20°C.

## CAUTION

The reagent contains phenol. Handle with care. The safety statements are on the label. It is advised to read SDS before reagent handling. Waste products must be handled as per local regulations.

## QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

## AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

## PROCEDURE

Procedure	BL mL	SA mL	ST mL
Sample	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min at 37°C, or 10 min at room temperature (20 - 25°C).

## Reading

Wavelength: 546 nm; 505 nm  
Blank: the contents of BL  
Colour stability: 1 hour (when protected from direct sunlight).

## CALCULATIONS

$$\frac{\text{SA Abs.}}{\text{ST Abs.}} \times 200 = \text{mg Cholesterol/dL}$$

Where:

SA Abs: Sample Absorbance  
ST Abs: Standard Absorbance

## SI Units

$$(\text{mg}/100 \text{ dL}) \times 0.0259 = \text{mmol/L}$$

## REFERENCE VALUES

### Coronary heart disease risk:

< 200 mg/dL	Desirable
200 - 239 mg/dL	Borderline high
> 239 mg/dL	High

## PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method.

Sensitivity, as detection limit: 2 mg/dL

Linearity: 700 mg/dL. For higher concentrations dilute the sample 1/2 with saline (NaCl 0.9%). Multiply the final result by 2.

Accuracy: 98.6%

Repetitivity, as CV%: 0.87%

Reproducibility, as CV%: 1.44%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request

## INTERFERENCES

There are no interferences by hemoglobin up to 200 mg/dL, as well as by bilirubin up to 15 mg/dL.

Do not pipette directly from the bottle of reagent, to avoid undesirable contamination.

## REFERENCES

- Fiedewald, W., Levy, R., Fredrickson, D.S. (1972), Clin.Chem, 18, 499-502.  
Allain, C.C., Poon, L.S., Chan, C.S.G., Richmond, W., Fu, P.C. (1974) Clin. Chem., 20, 470-475.  
Zoppi, F., Fellini, D. (1976), Clin. Chem., 22, 690-691.  
Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).  
CLSI Guidelines and Standards, CLSI, Wayne, PA  
Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).  
National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.  
Zapico E. y Ordóñez J. (2002), Clin Invest Ar ter ioscl;14(5):272-6.

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### BIBLIOGRAPHIE

- Fiedewald, W., Levy, R., Fredrickson, D.S. (1972), Clin.Chem, 18, 499-502.  
Allain, C.C., Poon, L.S., Chan, C.S.G., Richmond, W., Fu, P.C. (1974) Clin. Chem., 20, 470-475.  
Zoppi, F., Fellini, D. (1976), Clin. Chem., 22, 690-691.  
Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).  
CLSI Guidelines and Standards, CLSI, Wayne, P.A  
Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).  
National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication.  
Bethesda: National Heart, Lung, and Blood Institute.  
Zapico E. y Ordóñez J. (2002), Clin Invest Ar ter ioscl, 14(5):272-6.

---

# COLESTEROL LIQUIDO CHOLESTEROL LIQUID

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Contiene: Dodecano-1-ol, etoxilado

GB - CHOLESTEROL ENZYMATIC, STANDARD 200 MG/DL

### Danger

Hazard: Causes serious eye damage.

Precautionary: Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.

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