TRIGLICÉRIDOS LÍQUIDOS

MÉTODO GPO

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Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACCPress (2000).

National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

Guidelines for the management of dyslipidaemias. (2016). European Heart Journal, 37, 2999-3058.


TRIGLYCERIDES LIQUID

GPO METHOD

For "in vitro" determination of triglycerides in serum plasma

PRINCIPLE

Triglycerides present in the sample are enzymatically hydrolysed by the action of lipases leading to glycerol and fatty acids. In presence of glycerol kinase (GK), the ATP phosphorylates glycerol to give glycerol-3-phosphate and the corresponding ADP. By glycerophosphate oxidase (GPO), glycerol -3 - phosphate is oxidized to dihydroxyacetone phosphate and hydrogen peroxide. In the last stage, with the peroxidase as a catalyst, hydrogen peroxide reacts with 4-aminoantipyrine and 4-chlorophenol to give quinoneimine. The intensity of the red colour is proportional to the amount of triglycerides present in the sample.

DIAGNOSTIC USE

The increased blood triglyceride level is a risk factor in the development of coronary heart disease. About 50% of the lipids of atheromatous lesions occurring in the coronary arteries are triglycerides, making it possible to relate to the pathogenesis of coronary arteriosclerosis.

The determination of triglycerides allows assessing the early risk of developing coronary atherosclerosis. The TGL may be increased by chylomicronemia (Fredrickson types I and V), type IV hypertriglyceridemia, diabetes mellitus, insulin resistance, obesity, hypothyroidism, pancreatitis, glycogen storage disease, nephrotic syndrome, hypertriglyceridemia sensitive to carbohydrates, Tangier's disease, von Gierke's disease, pernicious anemia, acute pancreatitis, Down's syndrome, biliary cirrhosis, septicemia.

Single test result cannot be used to make a clinical diagnosis. The results are to be evaluated in the context of all clinical and laboratory data obtained.

REAGENTS

Kit 1 x 100 mL. (Ref. 99 23 30) Contents:

A. 1 x 100 mL Reagent

Ref. 99 23 25

B. 1 x 5 mL Standard

Ref. 99 03 17

Kit 3 x 100 mL. (Ref. 99 23 20). Contents:

A. 3 x 100 mL Reagent

Ref. 99 23 25

B. 1 x 5 mL Standard

Ref. 99 03 17

Kit 2 x 250 mL. (Ref. 99 30 80). Contents:

A. 2x 250 mL Reagent

Ref. 99 01 53

B. 1 x 5 mL Standard

Ref. 99 03 17

WORKING REAGENT PREPARATION

Reagent and standard are ready to use.

REAGENT COMPOSITION

Concentration in the reagent is:

Pipes buffer pH 6,8	50 mM
4-chlorophenol	4.2 mM
4-aminoantipyrine	0,35 mM
ATP	2 mM
Magnesium aspartate	40 mM
Glycerol kinase	≥ 800 U/l
Glycerol-3-phosphate oxidase	≥ 2000 U/l
Peroxidase	≥ 500 U/l
Lipases	≥ 9000 U/l
Non reactive stabilizers	

Standard: Solution of glycerol in water equivalent to 200 mg/dL (2.29 mmol/L)

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank >0.400

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C. Cuvette: 1 cm light-path

SAMPLE

Serum, or plasma with EDTA or heparin. Triglycerides are stable 4 days at 2-8°C, and up to three months at -20°C

CAUTION

The reagent contains phenol derivatives. Handle with care.

The safety statements are on the label. It is advisable to look at the SDS before using the reagent.

The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

It is recommended to include control sera, Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85), in each test series for assuring the obtained values.

Each particular laboratory should establish its own quality control program and measures for avoiding results deviations

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

Bring the reagent and the analyzer to working temperature.

Technique	BL	SA	ST
	mL	mL	mL
Standard	---	---	0.01
Sample	---	0.01	---
Reagent	1.00	1.00	1.00

Mix well and let stand 5 min at 37°C or 10 min. at room temperature (15 - 25°C)

Reading

Wavelength: 546 nm; 505 nm

Blank: the contents of BL

Colour stability: 1 hour, when protected from direct sunlight

CALCULATIONS

$$\frac{\text{SA Abs}}{\text{ST Abs}} \times 200 = \text{mg of triglycerides/dL}$$

Where:

SA Abs: Sample absorbance

ST Abs: Standard absorbance

S.I. Units

$$(\text{mg/dL}) \times 0.01143 = \text{mmol/L}$$

REFERENCES VALUES

According to the recommendations of the European Atherosclerosis Society, Europea Society of Cardiology and NCEP:

Normal: <150 mg/dL (<1,7 mmol/L)

Borderline: 150 - 199 mg/dL (1,7 - 2,25 mmol/L)

High: 200 - 499 mg/dL (2,26 - 5,64 mmol/L)

Very high: >500 mg/dL (>5,64 mmol/L)

PERFORMANCE CHARACTERISTICS

Performance of the reagent depends on the reagent itself and also depends on method and analyzer used. The results indicated are obtained using a manual method.

Sensitivity, as detection limit: 3.0 mg/dL

Linearity: Up to 1000 mg of Triglycerides/dL. Samples with a higher concentrations shall be diluted 1/10 with NaCl 0.9% and assayed once again. Multiply the final result by 10.

Accuracy: 98.5 %.

Repetitivity as Variation Coefficient: 0.89%

Reproducibility as Variation Coefficient: 1.52%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

INTERFERENCES

The assay is not interfered by hemoglobin (up to 150 mg/dL) nor bilirubin (up to 20mg/dL).

The use of disposable equipment is recommended to prevent unwanted contamination.

REFERENCES

Jacobs,N.J., VanDemark,P.J. (1960). J. Bacteriol.,79, 532 - 538.

Trinder,P. (1969). Ann.Clin. Biochem.,6, 24 - 27.

A policy statement of the European Atherosclerosis Society, European Heart Journal 8,(1987) 77 - 88.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Educacion Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

Guidelines for the management of dyslipidaemias. (2016). European Heart Journal, 37, 2999-3058.

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